


APPLICATION




OCT works in the mesoscale and bridges the gap between microscopic methods like Confocal Laser Scanning Microscopy (CLSM) and macroscopic methods like photography.¹

QUICK FACTS

- ◆ OCT is a non-destructive, non-invasive imaging technique.
- ◆ Biofilm imaging is usually performed in water. Because spectral absorption in water is lower at 880 nm than at 1300 nm, OCT systems operating near 880 nm are preferred.
- ◆ No dyes are needed; OCT works by analyzing backscattered light.
- ◆ Imaging through a glass or PMMA window (such as in a flow cell) is possible.
- ◆ Typical interests include morphology and changes in morphology over time due to growth or external forces.
- ◆ OCT can image into biofilms and aids in the analysis of structural features such as pore density and size.
- ◆ Flow analysis is possible via Doppler OCT (Doppler mode included in all systems).

TYPICAL SETUP

- ◆ The OCT scan lens can be dipped in open water for imaging.
- ◆ The biofilm can be grown in a reactor and then be taken out and analyzed under the OCT scanner.
- ◆ Imaging in a flow cell setup is possible, as illustrated to the right.




PUBLICATIONS



- 1) M. Wagner, D. Taherzadeh, C. Haisch, H. Horn, *Biotechnol. Bioeng.*, **107** (5), 844, 2010
- 2) L. Fortunato, S. Jeong, T.O. Leiknes, *Sci. Rep.*, **7**, 15, 2017
- 3) C. Picioreanu, F. Blauert, H. Horn, M. Wagner, *Water Res.*, **145**, 588, 2018
- 4) K.J. Martin, D. Bolster, N. Derlon, E. Morgenroth, R. Nerenberg, *J. Membr. Sci.*, **471**, 130, 2014
- 5) Y. Gao, S. Haavisto, C.Y. Tang, J. Salmela, W. Li, *J. Membr. Sci.*, **448**, 198, 2013
- 6) C. Dreszer, A.D. Wexler, S. Drusová, T. Overdijk, A. Zwijnenburg, H.C. Flemming, J.C. Kruithof, J.S. Vrouwenvelder, *Water Res.*, **67**, 243, 2014
- 7) N. Derlon, M. Peter-Varbanets, A. Scheidegger, W. Pronk, E. Morgenroth, *Water Res.*, **46**, 3323, 2012
- 8) M. Wagner, H. Horn, *Biotechnol. Bioeng.*, **114** (7), 1386, 2017
- 9) A.F. Rosenthal, J.S. Griffin, M. Wagner, A.I. Packman, O. Balogun, G.F. Wells, *Biotechnol. Bioeng.*, **115** (9), 2268, 2018

* Images licensed under a Creative Commons Attribution 4.0 International License, see <http://creativecommons.org/licenses/by/4.0/>. Image adapted from Original.


EXAMPLE IMAGES



Growth of biofilm on a membrane over several days. The accumulation of biofilm on a membrane leads to a significant reduction in performance and is therefore of interest in many industrial applications.^{2,*}




Deformation of Biofilm Due to Shear Forces: (a) Original Outline of the Biofilm, (b) Deformation Under Flow of Water from Left to Right, (c) Sketch Comparing the Structure With and Without Applied Force^{3,*}



RECOMMENDED ITEMS

Choice of OCT System:

- ◆ **GAN312C1(/M)**: For Standard Imaging of Thick Biofilms
- ◆ **GAN612C1(/M)**: For Imaging Dynamic Processes
- ◆ **GAN332C1(/M) or GAN632C1(/M)**: For High-Resolution Imaging of Thin Biofilms



Useful Accessories:

- ◆ For Larger Depth of Focus: **OCT-LK4-BB** Scan Lens Kit and **OCT-RA4** Reference Length Adapter
- ◆ Recommended When Imaging Through More Than 10 mm of Water, Glass, and PMMA: **SRA10** Spacer Reference Arm
- ◆ Special Spacers and Reference Arms for Imaging in Open Water:
 - **OCT-IMM3-SP1 & SRA10** for OCT-LK3-BB Lens
 - **OCT-IMM4-SP1 & SRA30** for OCT-LK4-BB Lens