


Brain Angiography

APPLICATION



(a) Maximum Intensity Projection (MIP) of blood vessel network in a mouse brain. (b) 3D red blood cell flux of the red square in (a) extracted via deep learning.^{1,*}

OCT Angiography uses the signal from moving blood cells to highlight blood vessels amongst the surrounding static tissue. Thus, it provides a useful tool for understanding brain damage caused by cerebrovascular dysfunctions.

QUICK FACTS



- ◆ OCT Angiography highlights blood vessels through changes in the OCT signal caused by moving blood cells.
- ◆ No dyes are necessary.
- ◆ OCT Angiography has to be performed *in vivo*.
- ◆ The Speckle Variance Angiography Mode is included in the complimentary ThorImage®OCT software package.
- ◆ Typical penetration depths are 1.5 to 2 mm in brain tissue.

TYPICAL SETUP

- ◆ Restraining the head reduces motion artifacts.
- ◆ It is recommended to thin the skull or use optical clearing to reduce signal losses due to absorption and scattering.²

- ◆ A cranial window preserves the integrity of the skull.^{3,4,5,6,7}

- ◆ A novel long-term clearing cranial window (LCCW) can preserve the skull transparency for 2 months with a single surgery.²



Representative photograph of dry skull (Initial), optical clearing window (Clearing), and LCCW in eight-week-old mice.^{2,*}

RECOMMENDED ITEMS

Choice of OCT System:

- ◆ High Resolution: **TEL221C1 (M)**
- ◆ High Penetration Depth: **VEG210C1 (M)**




TEL221 C1


Useful Accessories:

- ◆ Different Objectives for Different Purposes:
 - High-Resolution Objective **OCT-LK2** for Small Capillary Imaging
 - Long-Focus Objective **OCT-LK4** for Large Depth of Focus (Deep Imaging)

EXAMPLE IMAGES



(a) *En face* Video Image of Mouse Brain (b) OCT cross section image. (c-f) Mouse blood vessel images acquired using speckle variance OCT. Thick blood vessels are found at the top of the brain (c), and thinner capillaries are located deeper (e-f).^{**}



Imaging blood vessels through dry skull (Initial), optical clearing window (Clearing), and LCCW in eight-week-old mice. The maps are vascular maximum intensity projection (MIP) views in the depth direction (from surface to 560 μm depth). The cross-sectional angiograms are from yellow dotted lines (Gray: OCT structure; Red hot: blood vessels)^{2,*}

PUBLICATIONS

- 1) S. Stefan, A. Kim, P. J. Marchand, F. Lesage, J. Lee, *Front. Neurosci.*, **16**, 835773, 2022.
- 2) C. Zhang, C.-J. Liu, W. Feng, *Adv. Sci.*, **9** (17), 2105893, 2022.
- 3) J. Tang, S.E. Erdener, S. Sunil, D.A. Boas, *J. Biomed. Opt.*, **24** (3), 036005, 2019.
- 4) W.J. Choi, R.K. Wang, *J. Biomed. Opt.*, **20** (10), 106004, 2015.
- 5) A. Akif, K. Walek, C. Polucha, J. Lee, *Biomed. Opt. Express.*, **9** (11), 5340, 2018.
- 6) C.W. Merkle, V.J. Srinivasan, *NeuroImage*, **125**, 350, 2016.
- 7) Á. Nyúl-Tóth, S. Tarantini, J. Delfavero, F. Yan, P. Balasubramanian, A. Yabluchanskiy, C. Ahire, T. Kiss, T. Csipo, A. Lipecz, A. E. Farkas, I. Wilhelm, I. A. Krizbai, Q. Tang, A. Csizsar, Z. Ungvari, *Am. J. Physiol. Heart Circ. Physiol.* **320** (4), H1370, 2021.

* Images licensed under a Creative Commons Attribution 4.0 International License, see <http://creativecommons.org/licenses/by/4.0/>. Image adapted from original.

** Images Acquired in Collaboration with MacVicar Lab, University of British Columbia.

Interested? Email OCT@thorlabs.com for more information.