



Date	9/8/2021	Production Order Number	TP02112949
Serial Number	SN: 00001	Assembled By	MM
Range Tested	25mm	Technician	CS



Pitch Deviation ^a (µRad)	
Min	-38.50
Max	17.40
Range	55.90

^a 25 mm Test Range




Yaw Deviation ^a (µRad)	
Min	-0.80
Max	46.20
Range	47.00

^a 25 mm Test Range




Date	9/8/2021	Production Order Number	TP02112949
Serial Number	SN: 00001	Assembled By	MM
Range Tested	25mm	Technician	CS



Pitch Deviation ^a (µRad)	
Min	-58.80
Max	0.00
Range	58.80

^a. 25 mm Test Range




Yaw Deviation ^a (µRad)	
Min	-30.70
Max	3.90
Range	34.60

^a. 25 mm Test Range




Date	9/8/2021	Production Order Number	TP02112949
Serial Number	SN: 00001	Assembled By	MM
Range Tested	25mm	Technician	CS



Pitch Deviation ^a (µRad)	
Min	-37.50
Max	2.90
Range	40.40

^a. 25 mm Test Range



Yaw Deviation ^a (µRad)	
Min	-78.90
Max	0.00
Range	78.90

^a. 25 mm Test Range

