



# PDX1AV(/M) Test Report

Date	4/2/2024	Production Order Number	W123456
Serial Number	12345601	Assembled By	LQ_Jiang
Range Tested	20mm	Techician	LQ_Jiang


### Forward Pitch ( $\mu\text{rad}$ ) vs. Position (mm)




### Backward Pitch ( $\mu\text{rad}$ ) vs. Position (mm)



**Forward Yaw ( $\mu\text{rad}$ ) vs. Position (mm)**



**Backward Yaw ( $\mu\text{rad}$ ) vs. Position (mm)**



Items	Pitch/ $\mu\text{Rad}$	Yaw/ $\mu\text{Rad}$
Forward	70	85
Backward	67	88
Specification	200	200